

3076120 Clino® STE7s Cleaning trolley for Healthcare

**Individualization
on request!**



Product description

Open trolley with tubSieve method, disposal and surface cleaning.

- Cleaning method: TubSieveMethod
- Mop change method for healthcare cleaning
- Ready-to-use delivery and tool-free customization



Product specifications

Article name Article number	Clino® STE7s 3076120
Weight	25.06 kg
Working dimension	107 x 61 x 115 cm
Resting dimension	107 x 61 x 115 cm
Basic components	Module Clino® STE 8/9 4x Clino® plastic swivel castor 125 mm 4x deflector ring Bucket 25 L blue Bucket 15 L Ergo blue with handle ClinoBar bottle holder blue Impact protection Draining sieve Bowl Extend blue Clino® ST/STE ErgoGrip blue Bucket 6 L red Bucket 6 L blue Bucket 6 L yellow Bucket 6 L green Device holder 22 mm Clino® STE bag holder frame Insert x 120 L Clino® STE Clino® STE bag support Clino® ST/STE ErgoGrip blue
Work area	The cleaning trolley should be completely cleaned after use and stored in a dry place. The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Supply area	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Disposal area	1 Pieces 1 PU
Care instructions	1 Pieces
Storage	25.06 kg
Pieces per PackingUnit Selling unit	no pallet shipment 120 x 80 x 220 cm
Pieces per box	87168000
Carton weight	Germany
Cartons per pallet Pallet	4065237017975
Customs tariff number	Our sales department will be happy to advise you on suitable accessories
Country of origin	
EAN-Number	
Compatible accessories	

Industry Sectors Assembly



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



3076120 Clino® STE7s

Cleaning trolley for Healthcare

1x	3070420	Module Clino® STE 8/9
4x	2013115	Clino® plastic swivel castor 125 mm
4x	2014030	Deflector ring - for pipe 22 mm
1x	3070100	Clino® STE bag holder frame with hook
1x	3070120	Folding disposal unit 1 x 120 litres
1x	3070200	Clino® STE bag support
3x	2800002	Clino® ST/STE ErgoGrip blue
1x	2870003	Clino® Bar bottle holder blue
3x	2808012	Bowl Extend blue
1x	2017100	Drip sieve
1x	2025002	Bucket 25 L blue
1x	2815012	Bucket 15 L Ergo blue with handle
1x	2000601	Bucket 6 L red
1x	2000602	Bucket 6 L blue
1x	2000603	Bucket 6 L yellow
1x	2000604	Bucket 6 L green
1x	2017031	Device holder 22 mm
2x	3020630	Impact protection
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig.reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

