

3076270 Clino® STE7s Cleaning trolley for Healthcare

**Individualization
on request!**



Product description

Open trolley with tubSieve method, disposal and surface cleaning.

- Cleaning method: TubSieveMethod
- Mop change method for healthcare cleaning
- Ready-to-use delivery and tool-free customization



Product specifications

Article name Article number	Clino® STE7s 3076270
Weight	24.2 kg
Working dimension	155 x 61 x 115 cm
Resting dimension	107 x 61 x 115 cm
Basic components	Basic module (cleanroom-suitable stainless steel, tube thickness 22 mm) Double ball bearing swivel castors 125 mm ErgoGrip-R Stainless steel drip tray Clino® Bar bottle holder Bucket 25 L-R Bucket 15 L Ergo-R Bucket 6 LR Shells Extend-R Device holder Unit 2 x 70 L
Work area	
Supply area	
Disposal area	
Care instructions	The cleaning trolley should be completely cleaned after use and stored in a dry place. The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Storage	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton size Carton weight	110 x 60 x 118 cm 24.2 kg
Cartons per pallet Pallet	1 120 x 80 x 220 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237017982
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



3076270 Clino® STE7s

Cleaning trolley for Healthcare

1x	3070420	Module Clino® STE 8/9
4x	2013115	Clino® plastic swivel castor 125 mm
4x	2014030	Deflector ring - for pipe 22 mm
1x	3070100	Clino® STE bag holder frame with hook
1x	3070270	Folding disposal unit 2 x 70 litres
1x	3070200	Clino® STE bag support
3x	2800002	Clino® ST/STE ErgoGrip blue
1x	2870003	Clino® Bar bottle holder blue
3x	2808012	Bowl Extend blue
1x	2017100	Drip sieve
1x	2025002	Bucket 25 L blue
1x	2815012	Bucket 15 L Ergo blue with handle
1x	2000601	Bucket 6 L red
1x	2000602	Bucket 6 L blue
1x	2000603	Bucket 6 L yellow
1x	2000604	Bucket 6 L green
1x	2017031	Device holder 22 mm
2x	3020630	Impact protection
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig.reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

