

3076300-R Clino® STEC7s-R Cleaning trolley for Healthcare

**Individualization
on request!**



Product description

Cleaning trolley with tub sieve method, disposal and surface cleaning made of recycled material.

- Cleaning method: TubSieveMethod
- Mop change method for healthcare cleaning
- Plastic components made from up to 99% recycled material, stainless steel module made from 80% recycled raw materials
- Ready-to-use delivery and tool-free customization



Product specifications

Article name Article number	Clino® STEC7s-R 3076300-R
Weight	27.7 kg
Working dimension	155 x 61 x 115 cm
Resting dimension	107 x 61 x 115 cm
Basic components	Basic module (cleanroom-suitable stainless steel, tube thickness 22 mm) Double ball bearing swivel castors 125 mm ErgoGrip-R Stainless steel drip tray Clino® Bar bottle holder Bucket 25 L-R Bucket 15 L Ergo-R Clino® MultiBoxes-R Clino® PlusBoxes Device holder Unit 2 x 70 L
Work area	62,34 %
Supply area	The cleaning trolley should be completely cleaned after use and stored in a dry place. The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Disposal area	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Recycling rate	1 Pieces 1 PU
Care instructions	1 Pieces
Storage	110 x 60 x 118 cm 27.7 kg
Pieces per PackingUnit Selling unit	1 120 x 80 x 220 cm
Pieces per box	87168000
Carton size Carton weight	Germany
Cartons per pallet Pallet	4065237034378
Customs tariff number	Our sales department will be happy to advise you on suitable accessories
Country of origin	
EAN-Number	
Compatible accessories	

Industry Sectors Assembly Recycling



Healthcare



Facility



Ready-to-use



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



3076300-R Clino® STEC7s-R Cleaning trolley for Healthcare

1x	3070420	Module Clino® STE 8/9
4x	2013115	Clino® plastic swivel castor 125 mm
4x	2014030	Deflector ring - for pipe 22 mm
1x	3070100	Clino® STE bag holder frame with hook
1x	3070270	Folding disposal unit 2 x 70 litres
1x	3070200	Clino® STE bag support
3x	7800006	ErgoGrip-R anthracite
1x	7870006	Clino® Bar-R bottle cage anthracite
1x	7088665	Clino® MultiBox-R anthracite centre/bottom
2x	7088660	Clino® MultiBox-R anthracite top
1x	2088680	Guide rails for Clino® MultiBox centered
1x	2088685	Guide rails for Clino® MultiBox bottom
1x	2017100	Drip sieve
1x	7025002	Bucket 25 L-R black Handle blue
1x	7815012	Bucket 15 L Ergo-R anthracite Handle blue
1x	2088611	Clino® PlusBox red
1x	2088622	Clino® PlusBox blue
1x	2088633	Clino® PlusBox yellow
1x	2088644	Clino® PlusBox green
4x	2088650	Transparent lid
4x	2088653	Transparent use
1x	2088655	Sieve transparent
1x	2017031	Device holder 22 mm
2x	3020630	Impact protection
1x	3060936	Clino® key
1x	9990810	PPS lanyard anthracite
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

