

3076300 Clino® STEC7s Cleaning trolley for Healthcare

Individualization
on request!



Product description

Cleaning trolley with tub sieve method, disposal and surface cleaning.

- Cleaning method: TubSieveMethod
- Mop change method for healthcare cleaning
- Ready-to-use delivery and tool-free customization

Product specifications

Article name Article number	Clino® STEC7s 3076300
Weight	27.7 kg
Working dimension	155 x 61 x 115 cm
Resting dimension	107 x 61 x 115 cm
Basic components	Basic module (cleanroom-suitable stainless steel, tube thickness 22 mm) Double ball bearing swivel castors 125 mm ErgoGrip-R Stainless steel drip tray Clino® Bar bottle holder Bucket 25 L-R Bucket 15 L Ergo-R Clino® MultiBoxes-R Clino® PlusBoxes Device holder Unit 2 x 70 L
Work area	The cleaning trolley should be completely cleaned after use and stored in a dry place.
Supply area	The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Disposal area	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Care instructions	
Storage	
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton size Carton weight	110 x 60 x 118 cm 27.7 kg
Cartons per pallet Pallet	1 120 x 80 x 220 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237017999
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



3076300 Clino® STEC7s

Cleaning trolley for Healthcare

1x	3070420	Module Clino® STE 8/9
4x	2013115	Clino® plastic swivel castor 125 mm
4x	2014030	Deflector ring - for pipe 22 mm
1x	3070100	Clino® STE bag holder frame with hook
1x	3070270	Folding disposal unit 2 x 70 litres
1x	3070200	Clino® STE bag support
3x	2800002	Clino® ST/STE ErgoGrip blue
1x	2870003	Clino® Bar bottle holder blue
1x	2088625	Clino® MultiBox blue
2x	2088620	Clino® MultiBox blue
1x	2088680	Guide rails for Clino® MultiBox centered
1x	2088685	Guide rails for Clino® MultiBox bottom
1x	2017100	Drip sieve
1x	2025002	Bucket 25 L blue
1x	2815012	Bucket 15 L Ergo blue with handle
1x	2088611	Clino® PlusBox red
1x	2088622	Clino® PlusBox blue
1x	2088633	Clino® PlusBox yellow
1x	2088644	Clino® PlusBox green
4x	2088650	Transparent lid
4x	2088653	Transparent use
1x	2088655	Sieve transparent
1x	2017031	Device holder 22 mm
2x	3020630	Impact protection
1x	3060936	Clino® key
1x	9990810	PPS lanyard anthracite
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

