

TECHNICAL DATASHEET

3079120 Clino® STE9L

Cleaning trolley for Healthcare, Facility

Individualization
on request!



Product description

Open trolley with pre-preparation method, disposal and surface cleaning.

- Cleaning method: EasyMop® pre-preparation method
- Maximum efficiency and safety in a single-stage process
- Ready-to-use delivery and tool-free customization

Product specifications

Article name Article number	Clino® STE9L 3079120
Weight	26.63 kg
Working dimension	107 x 61 x 115 cm
Resting dimension	107 x 61 x 115 cm
Basic components	Module Clino® STE 8/9 4x Clino® plastic swivel castor 125 mm 4x deflector ring Flat wringer RINGO ST red Bucket 25 L red Bucket 15 L Ergo blue with handle ClinoBar bottle holder blue Impact protection Bowl Extend blue Clino® ST/STE ErgoGrip blue Bucket 6 L red Bucket 6 L blue Bucket 6 L yellow Bucket 6 L green Device holder 22 mm Clino® STE bag holder frame Insert x 120 L Clino® STE Clino® STE bag support Clino® ST/STE ErgoGrip blue
Work area	The cleaning trolley should be completely cleaned after use and stored in a dry place. The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Supply area	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Disposal area	
Care instructions	
Storage	
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton weight	26.63 kg
Cartons per pallet Pallet	no pallet shipment 120 x 80 x 220 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237018439
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors Assembly



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



3079120 Clino® STE9L

Cleaning trolley for Healthcare, Facility

1x	3070420	Module Clino® STE 8/9
4x	2013115	Clino® plastic swivel castor 125 mm
4x	2014030	Deflector ring - for pipe 22 mm
1x	3070100	Clino® STE bag holder frame with hook
1x	3070120	Folding disposal unit 1 x 120 litres
1x	3070200	Clino® STE bag support
3x	2800002	Clino® ST/STE ErgoGrip blue
1x	2870003	Clino® Bar bottle holder blue
3x	2808012	Bowl Extend blue
1x	2920001	Tray for EasyMop® system box blue
1x	2920010	Tray for EasyMop® system box blue
2x	2910000	Systembox EasyMop® blue
1x	2910010	Systembox EasyMop® red
1x	2000601	Bucket 6 L red
1x	2000602	Bucket 6 L blue
1x	2000603	Bucket 6 L yellow
1x	2000604	Bucket 6 L green
1x	2017031	Device holder 22 mm
2x	3020630	Impact protection
1x	2920030	Measuring jug PPS transparent
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

