

TECHNICAL DATASHEET

3504300 Clino® CR Mini MF-GMP Cleaning trolley for Cleanroom

Individualization
on request!

Product description



Compact cleaning trolley with MopFloat® GMP dosing and impregnation system for highly sensitive areas made of stainless steel.

- o Maximum efficiency and process reliability in a single-stage process
- o Made from the highest quality stainless steel
- o High-end surface quality through electropolishing
- o Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR Mini MF-GMP 3504300
Weight	18.79 kg
Working dimension	63.6 x 53 x 67.2 cm
Resting dimension	63.6 x 53 x 52 cm
Basic components	Basic module (cleanroom-compatible stainless steel) High-end finish Autoclavable special castors 100 mm 1 MopFloat® GMP with dosing ring set 250 ml 1 stainless steel container 25 L 1 MopTen stainless steel container
Work area	The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Care instructions	The trolley's autoclave suitability depends on its individual components.
Sterility	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Storage	1 Pieces 1 PU 1 Pieces 18.79 kg no pallet shipment 120 x 80 x 220 cm
Pieces per PackingUnit Selling unit	87168000
Pieces per box	Germany
Carton weight	4065237119792
Cartons per pallet Pallet	Our sales department will be happy to advise you on suitable accessories
Customs tariff number	
Country of origin	
EAN-Number	
Compatible accessories	

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3504300 Clino® CR Mini MF-GMP

Cleaning trolley for Cleanroom

1x	3530166	CR Mini module - facelift	
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
1x	3530180	Adapter for stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3500810	MopFloat® GMP	Autoclavable (121 °C for 20 minutes)
1x	3500832	MopFloat® insert 250 ml	Autoclavable (121 °C for 20 minutes)
1x	3500032	Stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3530187	Stem for MopTen and EasyMop® system box	Autoclavable (121 °C for 20 minutes)
1x	3500022	MopTen	Autoclavable (121 °C for 20 minutes)
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

