

TECHNICAL DATASHEET

3504500 Clino® CR Mini EM-GMP2 Cleaning trolley for Cleanroom

Individualization
on request!

Product description

Small cleaning trolley with EasyMop GMP® pre-preparation method for highly sensitive areas made of stainless steel.



- o Maximum efficiency and process reliability in a single-stage process
- o Made from the highest quality stainless steel
- o High-end surface quality through electropolishing
- o Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR Mini EM-GMP2 3504500
Weight	17.74 kg
Working dimension	48.4 x 51.3 x 44.6 cm
Resting dimension	48.4 x 51.3 x 44.6 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable special rollers 100 mm 2 EasyMop GMP® system boxes with Mop Jump System The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open. The product has to be prepared for cleanroom sustainability before use. The trolley's autoclave suitability depends on its individual components.
Work area	
Care instructions	
Sterility	
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton size Carton weight	69 x 59 x 52 cm 17.74 kg
Cartons per pallet	no pallet shipment
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237119853
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



TECHNICAL DATASHEET

3504500 Clino® CR Mini EM-GMP2 Cleaning trolley for Cleanroom

1x	3530166	CR Mini module - facelift	
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
1x	3530154	Adapter for CR4	Autoclavable (121 °C for 20 minutes)
2x	3500200	EasyMop® GMP system box	Autoclavable (121 °C for 20 minutes)
2x	3500210	MopJump®	Autoclavable (121 °C for 20 minutes)
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

