

TECHNICAL DATASHEET

3504510 Clino® CR Mini EM-GMP1 Cleaning trolley for Cleanroom

Individualization
on request!

Product description

Small cleaning trolley with EasyMop GMP® pre-preparation method and disposal for highly sensitive areas made of stainless steel.



- o Maximum efficiency and process reliability in a single-stage process
- o Made from the highest quality stainless steel
- o High-end surface quality through electropolishing
- o Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR Mini EM-GMP1 3504510
Weight	15.49 kg
Working dimension	49.3 x 51.3 x 44.6 cm
Resting dimension	49.3 x 51.3 x 44.6 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable special rollers 100 mm 1 EasyMop GMP® system box with Mop Jump System 1 MopTen stainless steel container for mop discharge
Work area	The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Care instructions	The product has to be prepared for cleanroom sustainability before use. The trolley's autoclave suitability depends on its individual components.
Sterility	1 Pieces 1 PU
Pieces per PackingUnit Selling unit	1 Pieces
Pieces per box	69 x 59 x 52 cm 15.49 kg
Carton size Carton weight	no pallet shipment
Cartons per pallet	87168000
Customs tariff number	Germany
Country of origin	4065237119884
EAN-Number	Our sales department will be happy to advise you on suitable accessories
Compatible accessories	

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3504510 Clino® CR Mini EM-GMP1 Cleaning trolley for Cleanroom

1x	3530166	CR Mini module - facelift	
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
1x	3500210	MopJump®	Autoclavable (121 °C for 20 minutes)
1x	3500022	MopTen	Autoclavable (121 °C for 20 minutes)
1x	3500200	EasyMop® GMP system box	Autoclavable (121 °C for 20 minutes)
1x	3530154	Adapter for CR4	Autoclavable (121 °C for 20 minutes)
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

