

TECHNICAL DATASHEET

3514100 Clino® CR1 FP-GMP Cleaning trolley for Cleanroom

Individualization
on request!



Product description

Cleaning trolley with Ringo® GMP flat wringer and disposal for highly sensitive areas made of stainless steel.

- Full flexibility in the one- and two-stage process
- Made from the highest quality stainless steel
- High-end surface quality through electropolishing
- Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR1 FP-GMP 3514100
Weight	17.4 kg
Working dimension	74 x 53 x 70 cm
Resting dimension	74 x 53 x 70 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable special rollers 100 mm Ringo® GMP flat wringer Stainless steel container 25 L Stainless steel container MopTen with Mop Jump System Device holder
Work area	Unit 1 x 35 L
Disposal area	The cleaning trolley should be completely cleaned after use and stored in a dry place.
Care instructions	The trolley's autoclave suitability depends on its individual components.
Sterility	1 Pieces 1 PU
Pieces per PackingUnit Selling unit	1 Pieces
Pieces per box	78 x 61 x 106 cm 17.4 kg
Carton size Carton weight	2 80 x 123 x 182 cm
Cartons per pallet Pallet	87168000
Customs tariff number	Germany
Country of origin	4065237108888
EAN-Number	Our sales department will be happy to advise you on suitable accessories
Compatible accessories	

Industry Sectors

Assembly



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3514100 Clino® CR1 FP-GMP Cleaning trolley for Cleanroom

1x	3530171	Module CR 1/2 - Facelift	Autoclavable (121 °C for 20 minutes)
1x	3530142	Bag support for CR1/CR4	Autoclavable (121 °C for 20 minutes)
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
1x	3530187	Stem for MopTen and EasyMop® system box	Autoclavable (121 °C for 20 minutes)
1x	3560000	Ringo® GMP flat wringer	Autoclavable (121 °C for 20 minutes)
1x	3530180	Adapter for stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3500032	Stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3500022	MopTen	Autoclavable (121 °C for 20 minutes)
1x	3500210	MopJump®	Autoclavable (121 °C for 20 minutes)
1x	2017031	Device holder 22 mm	
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

