

TECHNICAL DATASHEET

3514500 Clino® CR1 EM-GMP Cleaning trolley for Cleanroom

Individualization
on request!



Product description

Cleaning trolley with EasyMop GMP® pre-preparation method and disposal for highly sensitive areas made of stainless steel.

- o Maximum efficiency and process reliability in a single-stage process
- o Made from the highest quality stainless steel
- o High-end surface quality through electropolishing
- o Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR1 EM-GMP 3514500
Weight	15.6 kg
Working dimension	61 x 53 x 70 cm
Resting dimension	61 x 53 x 70 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable special rollers 100 mm EasyMop GMP® system boxes with Mop Jump System Device holder
Work area	Unit 1 x 35 L
Disposal area	The cleaning trolley should be completely cleaned after use and stored in a dry place.
Care instructions	The trolley's autoclave suitability depends on its individual components.
Sterility	
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton size Carton weight	78 x 61 x 106 cm 15.6 kg
Cartons per pallet Pallet	2 80 x 123 x 182 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237109007
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3514500 Clino® CR1 EM-GMP Cleaning trolley for Cleanroom

1x	3530171	Module CR 1/2 - Facelift	Autoclavable (121 °C for 20 minutes)
1x	3530142	Bag support for CR1/CR4	Autoclavable (121 °C for 20 minutes)
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
1x	3530154	Adapter for CR4	Autoclavable (121 °C for 20 minutes)
2x	3500200	EasyMop® GMP system box	Autoclavable (121 °C for 20 minutes)
2x	3500210	MopJump®	Autoclavable (121 °C for 20 minutes)
1x	2017031	Device holder 22 mm	
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

