

TECHNICAL DATASHEET

3514800 Clino® CR1 SC-CR Cleaning trolley for Cleanroom

Individualization
on request!



Product description

Cleaning trolley with MopScoop® CR dosing and impregnation system and disposal for non-sterile areas made of stainless steel and plastic.

- No downtime problems thanks to just-in-time impregnation
- Made from the highest quality stainless steel
- High-end surface quality through electropolishing
- Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR1 SC-CR 3514800
Weight	13.5 kg
Working dimension	78 x 53 x 70 cm
Resting dimension	78 x 53 x 70 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Electrically conductive special rollers 100 mm MopScoop® dosing unit incl. plastic lid EasyMop® system box Device holder
Work area	Unit 1 x 35 L
Disposal area	The cleaning trolley should be completely cleaned after use and stored in a dry place.
Care instructions	The trolley's autoclave suitability depends on its individual components.
Sterility	
Pieces per PackingUnit Selling unit	1 Pieces 1 PU
Pieces per box	1 Pieces
Carton size Carton weight	78 x 61 x 106 cm 13.5 kg
Cartons per pallet Pallet	2 80 x 123 x 182 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237109090
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3514800 Clino® CR1 SC-CR Cleaning trolley for Cleanroom

1x	3530171	Module CR 1/2 - Facelift	Autoclavable (121 °C for 20 minutes)
1x	3530142	Bag support for CR1/CR4	Autoclavable (121 °C for 20 minutes)
4x	2013430	Swivel castor 100 mm	
1x	3530187	Stem for MopTen and EasyMop® system box	Autoclavable (121 °C for 20 minutes)
1x	2910000	Systembox EasyMop® blue	
1x	2017031	Device holder 22 mm	
1x	3500710	MopScoop® GMP container	Autoclavable (121 °C for 20 minutes)
1x	3500722	Holding frame for MopScoop® GMP	
1x	3500730	Tray for MopScoop® GMP	Autoclavable (121 °C for 20 minutes)
2x	3500735	Storage tray for MopScoop® GMP tray	
1x	3500740	Trigger for MopScoop® GMP	Autoclavable (121 °C for 20 minutes)
1x	3500750	MopScoop® GMP dosing tube	Autoclavable (121 °C for 20 minutes)
1x	3500762	Nut for MopScoop® GMP	
1x	3500753	MopScoop® GMP dosing flask 200 ml	
2x	3500755	Seal for MopScoop® GMP dosing piston	
1x	3500756	Seal for MopScoop® GMP dosing tube	
1x	3500780	Lid for MopScoop® CR	
1x	0899120	Work instruction	
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

