

TECHNICAL DATASHEET

3544300 Clino® CR4 MF-GMP Cleaning trolley for Cleanroom

Individualization
on request!

Product description

Cleaning trolley with MopFloat® GMP and disposal for highly sensitive areas made of cleanroom-compatible stainless steel.

- Maximum efficiency and process reliability in a single-stage process
- Made from the highest quality stainless steel
- High-end surface quality through electropolishing
- Touch-free operating concept



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CR4 MF-GMP 3544300
Weight	26.4 kg
Working dimension	75 x 63 x 102 cm
Resting dimension	75 x 63 x 102 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable special rollers 100 mm with deflector ring MopFloat® GMP with dosing ring set 250 ml Stainless steel container 25 L MopTen stainless steel container with MopJump system Stainless steel container 4 L Handle holder Mop frame support
Work area	Unit 1 x 35 L The cleaning trolley should be completely cleaned after use and stored in a dry place. The trolley's autoclave suitability depends on its individual components.
Disposal area	1 Pieces 1 PU
Care instructions	1 Pieces
Sterility	78 x 61 x 106 cm 26.4 kg
Pieces per PackingUnit Selling unit	2 80 x 123 x 182 cm
Pieces per box	87168000
Carton size Carton weight	Germany
Cartons per pallet Pallet	4065237106785
Customs tariff number	Our sales department will be happy to advise you on suitable accessories
Country of origin	
EAN-Number	
Compatible accessories	

Industry Sectors

Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



3544300 Clino® CR4 MF-GMP Cleaning trolley for Cleanroom

1x	3530176	Module CR4 - Facelift	Autoclavable (121 °C for 20 minutes)
1x	3530184	Insert for bag holder frame	Autoclavable (121 °C for 20 minutes)
1x	3530142	Bag support for CR1/CR4	Autoclavable (121 °C for 20 minutes)
4x	3500100	Stainless steel swivel castor 100 mm	Autoclavable (121 °C for 20 minutes)
4x	3500084	Stainless steel cover plate	Autoclavable (121 °C for 20 minutes)
1x	3530180	Adapter for stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3500032	Stainless steel container 25 L	Autoclavable (121 °C for 20 minutes)
1x	3530187	Stem for MopTen and EasyMop® system box	Autoclavable (121 °C for 20 minutes)
1x	3500022	MopTen	Autoclavable (121 °C for 20 minutes)
1x	3500210	MopJump®	Autoclavable (121 °C for 20 minutes)
1x	3500093	Holder for handle	Autoclavable (121 °C for 20 minutes)
1x	3500092	Support for Mop frames	Autoclavable (121 °C for 20 minutes)
2x	3500094	Holding frame for stainless steel container 4 L	Autoclavable (121 °C for 20 minutes)
2x	3500007	Stainless steel container 4 L	Autoclavable (121 °C for 20 minutes)
1x	3500810	MopFloat® GMP	Autoclavable (121 °C for 20 minutes)
1x	3500832	MopFloat® insert 250 ml	Autoclavable (121 °C for 20 minutes)
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

