

TECHNICAL DATASHEET

3547000 Clino® CRX4 S-CR Cleaning trolley for Cleanroom

Individualization
on request!

Product description



Cleaning trolley with sieve, buckets and disposal for non-sterile areas made of stainless steel and plastic.

- Maximum efficiency and process reliability in a single-stage process
- Made of cleanroom-suitable stainless steel and plastic
- Touch-free operating concept
- This item is not supplied pre-assembled. Assembly instructions for self-assembly are included.



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® CRX4 S-CR 3547000
Weight	14.64 kg
Working dimension	75.7 x 60.4 x 100.4 cm
Resting dimension	75.7 x 60.4 x 100.4 cm
Basic components	Basic module (stainless steel) Tray for 3 15L buckets Electrically conductive special castors 100 mm
Work area	3 buckets 15 L 1 draining sieve 1 device holder
Disposal area	Unit 1 x 35 L
Care instructions	The cleaning trolley should be completely cleaned after use and stored in a dry place. The product has to be prepared for cleanroom sustainability before use. The trolley's autoclave suitability depends on its individual components.
Sterility	1 Pieces 1 PU
Pieces per PackingUnit Selling unit	1 Pieces
Carton size Carton weight	78 x 58 x 45 cm 14.64 kg
Cartons per pallet Pallet	8 80 x 123 x 182 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237120217
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;



TECHNICAL DATASHEET

3547000 Clino® CRX4 S-CR Cleaning trolley for Cleanroom

1x	3500231	CRX3 sheet metal	Autoclavable (121 °C for 20 minutes)
2x	3500240	Truss-head bolt CRX	Autoclavable (121 °C for 20 minutes)
4x	3500235	Spacer CRX	Autoclavable (121 °C for 20 minutes)
1x	3530184	Insert for bag holder frame	Autoclavable (121 °C for 20 minutes)
6x	3500239	Pan head screw CRX	Autoclavable (121 °C for 20 minutes)
2x	3500238	Stud screw CRX	Autoclavable (121 °C for 20 minutes)
3x	3550010	Bucket 15 L blue	Autoclavable (121 °C for 20 minutes)
1x	2017100	Drip sieve	
4x	3500236	CRX sealing washer	Autoclavable (121 °C for 20 minutes)
1x	3500093	Holder for handle	Autoclavable (121 °C for 20 minutes)
1x	3500233	Handle CRX4	Autoclavable (121 °C for 20 minutes)
4x	2013430	Swivel castor 100 mm	
1x	3500237	Allen key CRX 5 mm	Autoclavable (121 °C for 20 minutes)
1x	0899100	Cardboard small	
1x	0899911	Assembly instructions Clino® CRX	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

