

TECHNICAL DATASHEET

3547200 Clino® CRX4 FP-CR Cleaning trolley for Cleanroom

Individualization
on request!



Product description

Cleaning trolley with flat wringer and disposal made of stainless steel and plastic for non-sterile areas.

- Full flexibility in one- and two-stage processes
- Made from brushed stainless steel and plastic
- Touch-free operating concept
- This item is not supplied pre-assembled. Assembly instructions for self-assembly are included.



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name | Article number

Clino® CRX4 FP-CR | 3547200

Dimensions | Weight

76 x 60.5 x 100.5 cm | 17.8 kg

Working dimension

75.7 x 60.4 x 100.4 cm

Resting dimension

75.7 x 60.4 x 100.4 cm

Basic components

Basic module (stainless steel)

Tray for 3 buckets 15 litre

Electrically conductive special castors 100 mm

Work area

3 buckets 15 L

1 flat wringer

1 device holder

Disposal area

Unit 1 x 35 L

Care instructions

The cleaning trolley should be completely cleaned after use and stored in a dry place.

The product has to be prepared for cleanroom sustainability before use.

The trolley's autoclave suitability depends on its individual components.

Sterility

1 Pieces | 1 PU

Pieces per PackingUnit | Selling unit

1 Pieces

Pieces per box

Carton size | Carton weight

78 x 58 x 45 cm | 17.8 kg

Cartons per pallet | Pallet

8 | 120 x 80 x 220 cm

Customs tariff number

87168000

Country of origin

Germany

EAN-Number

4065237120279

Compatible accessories

Our sales department will be happy to advise you on suitable accessories

Industry Sectors



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

3547200 Clino® CRX4 FP-CR Cleaning trolley for Cleanroom

1x	3500233	Handle CRX4	Autoclavable (121 °C for 20 minutes)
1x	3530184	Insert for bag holder frame	Autoclavable (121 °C for 20 minutes)
1x	3500231	CRX3 sheet metal	Autoclavable (121 °C for 20 minutes)
2x	3500238	Stud screw CRX	Autoclavable (121 °C for 20 minutes)
2x	3500240	Truss-head bolt CRX	Autoclavable (121 °C for 20 minutes)
4x	3500235	Spacer CRX	Autoclavable (121 °C for 20 minutes)
4x	3500236	CRX sealing washer	Autoclavable (121 °C for 20 minutes)
6x	3500239	Pan head screw CRX	Autoclavable (121 °C for 20 minutes)
4x	2013430	Swivel castor 100 mm	
1x	2000070	Ringo® GMP Ultra flat wringer	Autoclavable (121 °C for 20 minutes)
3x	3550010	Bucket 15 L blue	Autoclavable (121 °C for 20 minutes)
1x	3500093	Holder for handle	Autoclavable (121 °C for 20 minutes)
1x	3500237	Allen key CRX 5 mm	Autoclavable (121 °C for 20 minutes)
1x	0899911	Assembly instructions Clino® CRX	
1x	0899100	Cardboard small	

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 30.01.2026;

