

7949270 Caro XT9-R

Cleaning trolley for Healthcare, Facility

**Individualization
on request!**



**MADE IN
EUROPE**

Product description

Cleaning trolley with pre-preparation method, storage space, disposal and surface cleaning made of recycled material.

- Cleaning method: EasyMop pre-preparation method
- Maximum efficiency and safety in a single-stage process
- Plastic components made from recycled material
- Ready-to-use delivery

Product specifications

Article name Article number	Caro XT9-R 7949270
Weight	33 kg
Working dimension	109 x 67 x 121 cm
Resting dimension	75 x 67 x 121 cm
Basic components	Basic module (plastic made from recycled material) Double ball bearing swivel castors 125 mm System boxes-R EasyMop® Half trays for Systembox-R EasyMop® Bucket XT-R 5 litre, handles Trays XT-R 20 L, extendable Tray XT-R 40 L, extendable Device holder Double hook
Work area	Unit XT-R 2 x 70 L
Supply area	Lid for disposal with SoftCloseSystem
Disposal area	78,16 % 5,052 kg
Recycling rate CO ² -savings	The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Care instructions	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Storage	1 Pieces 1 PU
Pieces per PackingUnit Selling unit	1 Pieces
Pieces per box	110 x 60 x 116 cm 33 kg
Carton size Carton weight	1 120 x 80 x 220 cm
Cartons per pallet Pallet	87168000
Customs tariff number	Germany
Country of origin	4065237035337
EAN-Number	Our sales department will be happy to advise you on suitable accessories
Compatible accessories	

Industry Sectors Assembly Hygiene Recycling



Healthcare



Facility



Ready-to-use



Touchless
operation



Components
up to
99%
recycled material

This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 05.01.2026;



7949270 Caro XT9-R

Cleaning trolley for Healthcare, Facility

1x	7990028	Intermediate shelf XT-R black
4x	5990026	Caro XT mounting bar
4x	5990045	Spacer sleeve black Caro XT
4x	5990035	Deflector ring grey Caro XT
2x	7920011	Tray for Systembox-R EasyMop®
4x	7990027	Side panel kit XT-R black
3x	7990052	Slot cover XT-R left black
3x	7990051	Slot cover XT-R right black
4x	2013114	Clino® plastic swivel castor 125 mm
1x	7990055	Bag holding frame XT-R without handle
2x	7990016	Double hook XT-R
1x	7990011	Bucket XT-R 7 L black, handle red
1x	7990012	Bucket XT-R 7 L black, handle blue
1x	7990013	Bucket XT-R 7 L black, handle yellow
2x	7990017	Tray XT-R 20 L extendable
1x	7990018	Tray XT-R 40 L extendable
1x	7990046	Device holder XT-R
2x	7990053	Slot cover XT-R basic module black
1x	7990021	Base plate XT-R black
2x	7990041	Plastic sleeve XT-R
1x	5990023	XT cover for frame
1x	7990034	Bag holding frame XT-R with handle
1x	5990036	Wheel disposal Caro XT
1x	7990022	Bag support XT-R black
1x	0899110	Carton large
1x	7910010	Systembox-R EasyMop® anthracite red
1x	7910000	Systembox-R EasyMop® anthracite blue
1x	5990004	Caro XT adapter for Clino® system boxes

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 05.01.2026;

