

7994270 Caro XT4-R

Cleaning trolley for Healthcare, Facility

**Individualization
on request!**



**MADE IN
EUROPE**

Product description

Cleaning trolley with EasyMop® pre-preparation method, disposal and surface cleaning from recycled material.

- Cleaning method: EasyMop® pre-preparation method
- Plastic components made from recycled material
- Maximum efficiency and safety in a single-stage process
- Ready-to-use delivery

Product specifications

Article name Article number	Caro XT4-R 7994270
Weight	23.2 kg
Working dimension	85 x 67 x 130 cm
Resting dimension	85 x 67 x 130 cm
Basic components	Basic module (plastic made from recycled material) Double ball bearing swivel castors 125 mm Half trays for Systembox-R EasyMop® Clino® system boxes-R Hanging basket with buckets 6 L-R Device holder Double hook
Work area	Unit XT-R 2 x 70 L Lid for disposal with SoftCloseSystem
Disposal area	76,52 %
Recycling rate	The cleaning trolley should be completely cleaned after use and stored in a dry place. The cleaning trolley should be completely cleaned after use and stored in a dry place. Store EasyMop system boxes open.
Care instructions	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Storage	1 Pieces 1 PU
Pieces per PackingUnit Selling unit	1 Pieces
Pieces per box	110 x 60 x 118 cm 23.2 kg
Carton size Carton weight	1 120 x 80 x 220 cm
Cartons per pallet Pallet	87168000
Customs tariff number	Germany
Country of origin	4065237035191
EAN-Number	Our sales department will be happy to advise you on suitable accessories
Compatible accessories	

Industry Sectors Assembly Hygiene Recycling



Healthcare



Facility



Ready-to-use



Touchless
operation



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;



7994270 Caro XT4-R

Cleaning trolley for Healthcare, Facility

1x	7990021	Base plate XT-R black
4x	7990053	Slot cover XT-R basic module black
4x	2013114	Clino® plastic swivel castor 125 mm
4x	5990045	Spacer sleeve black Caro XT
4x	5990035	Deflector ring grey Caro XT
2x	5990026	Caro XT mounting bar
2x	7990056	Side panel XT-R rounded black
3x	7990052	Slot cover XT-R left black
3x	7990051	Slot cover XT-R right black
1x	7990055	Bag holding frame XT-R without handle
1x	7990044	Insert 2 x 70 L
1x	7990048	Handle kit XT-R black
1x	5990023	XT cover for frame
1x	7990046	Device holder XT-R
2x	7990016	Double hook XT-R
1x	5990002	XT mounting frame for attaching accessories
1x	2019000	Hanging basket for three PlusBoxes or bucket 6L
1x	7000601	Bucket 6 L-R black Handle red
1x	7000602	Bucket 6 L-R black Handle blue
1x	7000603	Bucket 6 L-R black Handle yellow
1x	5990004	Caro XT adapter for Clino® system boxes
2x	7990041	Plastic sleeve XT-R
2x	7920011	Tray for Systembox-R EasyMop®
1x	7910000	Systembox-R EasyMop® anthracite blue
1x	7910010	Systembox-R EasyMop® anthracite red
1x	0899110	Carton large

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
 www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.11.2025;

